

Infinity Blanket FAQs

What is The Infinity Blanket made of?

The blanket is a blend of 72% Organic Cotton and 28% Celliant® Polyester.

Where is The Infinity Blanket made?

The Infinity Blanket is made in Germany from both organic and recycled materials.

What is Celliant®?

Celliant® is a revolutionary, patented fiber technology that utilizes a blend of natural minerals embedded in the fiber's core which convert the body's natural energy (heat) into infrared (IR) light which is emitted back into the tissue and muscles. This IR energy improves localized blood flow and increases the amount of oxygen reaching your muscles. This helps restore your muscles helping you recover faster, increases comfort and promotes restful sleep.

The FDA has determined that Celliant® products are medical devices as defined in section 201(h) of the Federal Food, Drug and Cosmetic Act and are general wellness products.

What is Celliant® made from?

Celliant® is a unique blend of 13 thermo-reactive minerals including titanium dioxide, silicon dioxide and aluminum oxide. Additional proprietary ingredients are blended with polyester fiber to create a variety of staple fibers, spun and filament yarns, and fabric blends.

Within each Infinity Blanket, approximately 12 plastic water bottles are recycled in the making of the Celliant® fibers.

Will Celliant®/SeaCell™ lose it's effectiveness over time?

No, the Celliant® will hold its strength over time (i.e., it will not wash out).

Does my Infinity Blanket have to touch my skin for a heightened effect?

Products enhanced with Celliant® technology can be covered or layered without losing effectiveness. Clothing, sheets, and laminated or waterproof fabrics do not hinder the efficacy of Celliant® to recycle energy to the body. The body's electromagnetic emissions travel through clothing, wood, glass and even concrete. The energy your body emits is the same that's seen in technology used by the military to view humans inside buildings.

Will I feel my Infinity Blanket working?

While you will not necessarily be able to perceive or feel an increase in circulation, oxygenation or blood flow, clinical testing has measured an increase in skin oxygen levels (TCPO2) and consistently recorded significant improvements across a wide spectrum of test subjects. Consumers who try Celliant® often experience increased comfort and an overall improved sense of health.

How do I care for my Infinity Blanket?

Wash cold and tumble dry low.

What is the size and weight of the Infinity Blanket?

The travel sized Infinity blanket is 105 x 150 cm big and weights 0,45 kg.

The throw sized Infinity blanket is 150 x 200 cm big and weights 0,9 kg.

Does the Infinity Blanket carry a certificate?

The organic cotton used for the Infinity Blanket is certified according to GOTS (Global Organic Textile Standard).